

Crosstex™ Patient Towels/Bibs/Liners

REF WPXAQ / WPXDR / WPXPE / WPXYE / WGXB L / W3D1930BL / WERUECBL
 WERUECGR / WPXAAQ / WPXADR / WPXAPE / WPXAYE / WPXBG / WPXGR
 WPXSI / WPXFP / WGXXWH / WERUECDR / WERUEC / WPXABG / WPXAGR
 WPXASI / WPXBL / WPXLV / WPXWH / WPXTB / WPXABL / WPXALV / WPXAWH
 WEXPP / WEXAQ / WEXDR / WEXPE / WEXYE / WUXBL / WEXAAQ / WEXADR
 WEXAPE / WEXAYE / WEXBG / WEXGR / WEXSI / WEXFP / WUXWH / WEXABG
 WEXAGR / WEXASI / WEXBL / WEXLV / WEXWH / WEXTB / WUXYE / WEXABL
 WEXALV / WEXAWH / WEPBAQ / WEPBGR / WEPBSI / WEPBTB / WTXAQ / WTXDR
 WTXYE / WT5BL / WEPBBL / WEPBLV / WEPBWH / WTXBG / WTXGR / WTXFP
 WT5WH / WEPBDR / WEPBPE / WEPBYE / WTXSI / WTXBL / WTXWH


Sani-Tab™ Chain-Free Patient Towels

REF WEXWHT / WEXLVT / WEXDRT / WEXBLT / WPXWHT / WPXDRT / WPXBLT

INSTRUCTIONS FOR USE

Crosstex™ Patient Towels/Bibs/Liners One towel per patient.

1. Attach clean clip to one side of towel.
2. Wrap around patient and attach other clip to towel.
3. Creating a barrier in front of the patient.



Sani-Tab™ Chain-Free Patient Towels One towel per patient.

1. Unfold towel and remove paper covering the two adhesive tabs from backside of towel.
2. Press adhesive tabs firmly against patient's collarbone to secure.
3. Creating a barrier in front of the patient.



CROSSTEX



Crosstex International, Inc.
 6789 W. Henrietta Road, Rush, NY 14543 USA
 585.359.0130 | 631.582.6777 | HuFriedyGroup.com



MDSS GmbH
 Schiffgraben 41
 30175 Hannover, Germany



UK REP STERIS Solutions Limited
 Chancery House, Rayns Way
 Watermead Business Park, Syston
 Leicestershire
 LE17 1PF
 United Kingdom

CH REP STERIS GmbH c/o BDO AG,
 Längfeldweg 116A,
 CH-2504 Biel-Bienne,
 SWITZERLAND



Made in USA

BIBS-IFU-EN / REV A - 12/2022

In case of any serious incident that has occurred and involves our devices within the European Union, please report it to Hu-Friedy Mfg. Co., LLC and your national competent authority of the member state.

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
 ©2022 Hu-Friedy Mfg. Co., LLC. All rights reserved.

Not made with natural rubber latex.